

PCR: When it doesn't work, add more stuff!

Reagent	Concentration	Purpose	Why it works	Caveats	Downstream issues
BSA	Up to 0.8 mg/ml	MAGIC (no seriously, just add some and things just work).	MAGIC (coats inside walls of tubes which prevents DNA from binding there, and also takes up space in tubes which can increase Taq performance.)		Too much can interfere with sequencing, but usually not.
Betaine	1-1.7M	Reduces secondary structures			
Formamide	1-5%	GC-rich regions, reduces secondary structures		Yes, it will cross-link the DNA a bit	Can inhibit Taq- you can add a touch more taq to compensate
DMSO	2-10% (start with 2%)	GC-rich regions, reduces secondary structures		10% will start killing off Taq	
Triton X-100	0.1-1%	Reduces secondary structures	Will increase yield	Increases non-specific amplification	
Tween-20/-40	9.5%	Reduces secondary structures Neutralizes SDS carryover		Increases non-specific amplifications	
MgCl₂	Already in master mix Can increase (or decrease) 1.0-4.0mM	Higher salt decreases Taq specificity. Lower salt increases specificity but eventually Taq won't function	Chelaters (EDTA, citrate) take up MgCl ₂ ; you might need to add more	Thaw and vortex- forms gradients in tubes over time.	Too much salt will carry over and inhibit subsequent sequencing reactions

TMAC	15-100mM	Increases primer specificity and melting temperature	Eliminates non-specific binding, helps with RNA contamination	Great for degenerate primers	Too much can interfere with sequencing reactions.
7-deaza-2'-deoxyguanosine (7-deaza dGTP)	Replace all G's with 7-deaza-2'-deoxyguanosine	Combats GC rich regions. Reduces duplex stability.		Can use 7-deanza dGTP:dGTP at 3:1 ratio.	Must be completely removed prior to sequencing reactions (sephadex or column purification recommended)